



# CLEANER DISINFECTANT Concentrated

# Jongro Podestoral Cleaner Disinfectant Cover tests Black 1 Model 6

### USE

Sanitiser and disinfectant.

### **DESCRIPTION**

Liquid disinfectant, deodoriser and multi-purpose cleaner. Kills a range of bacteria including MRSA and Shigella (cause of the most common type of Dysentery). Tested to kill E-Coli. For use in schools, hospitals, surgeries, nursing homes, sports clubs, swimming pools and wherever surfaces may be contaminated with harmful bacteria. Not suitable for food contact surfaces.

### **FEATURES & BENEFITS**

- Passes EN 1276 under dirty conditions at 1:25 with a 30 second contact time
- Kills MRSA and Shigella
- Helps prevent the spread of infection
- Leaves a clean fresh odour
- Non bleach formula
- Ideal in washrooms
- Suitable for use on a variety of hard surfaces

### **HOW TO USE**

### Disinfecting hard surfaces and equipment

Dilute 200ml per 5 litres of water (1:25) for mopping or wiping. For spray cleaning dilute 24ml per 600ml of water in a trigger bottle (1:25).

Allow 30 seconds contact time, then wipe away and rinse thoroughly.

# **General Cleaning**

Dilute 100ml per 5 litres of warm water (1:50) for mopping or wiping. For spray cleaning dilute 12ml per 600ml of water in a trigger bottle (1:50).

For heavy soiled surfaces repeat application.

# TYPICAL PRODUCT DATA

Appearance: Liquid BA048-5 2 x 5 Litre

Colour: Blue
Odour: Perfume
pH – undiluted: 11-12.5

# **BIODEGRADABILITY**

All surfactants used in this product comply with the current European Regulations concerning biodegradability and protection of the environment.

**PACKAGING**