

No: WPS/H/001

ISSUE: 19

DATE: Jan 2016

## 'hydrosoft' TABLET SALT

'hydrosoft' tablets are premium quality, additive free, compacted salt, one of INEOS Enterprises' range of superior products for water softener regeneration. The specification is assured by a Quality Management System registered to ISO 9001, so ensuring consistent performance.

## CHEMICAL SPECIFICATION

Test methods used are BS7319:1990 or equivalent, except appearance which is a visual assessment.

COMPONENT	UNIT	SPECIFICATION	TYPICAL ANALYSIS
Appearance		White	
Assay (dry basis)	%m/m NaCl	99.9 min.	99.9
Surface Moisture	%m/m H₂O	0.05 max.	0.02
Insoluble matter	mg/kg	<100	<10
*Alkalinity	mg/kgNa₂CO₃	<150	52
*Sulphate	mg/kgNa₂SO₄	<500	211
Total Iron	mg/kg Fe	<3	0.3
*Total Calcium	mg/kg Ca	<20	2.5
*Total Magnesium	mg/kg Mg	<5	0.6
*Total Copper	mg/kg Cu	2 max.	<0.1
*Total Arsenic	mg/kg As	0.3 max.	<0.01
*Total Lead	mg/kg Pb	1 max.	<0.1
*Total Cadmium	mg/kg Cd	0.2 max.	<0.01
*Total Mercury	mg/kg Hg	0.05 max.	<0.03
*Total Nickel	mg/kg Ni	0.75 max.	<0.05
*Total Chromium	mg/kg Cr	0.75 max.	<0.03
*Total Selenium	mg/kg Se	2.6 max.	<0.2
*Total Antimony	mg/kg Sb	2.6 max.	<0.2

<sup>\*</sup>measured on PDV

## PHYSICAL CHARACTERISTICS

Typical Tablet dimensions
Typical % m/m \\( \frac{1}{2} \)

Typical % m/m \\( \frac{1}{2} \)

Typical % m/m \( \frac{1}{2} \) 24 x 24 x 13 mm Approx.

Typical % m/m Whole Tablets > 80 Typical % m/m Fines < 2

> Important Note: The information contained in this document is given in good faith and is to the best of Darlingtons' knowledge correct at the date of publication, but it is for the users to satisfy themselves of the suitability of the product for their purposes.